

BS3+U

Handheld Dermatoscope Skin Analyzer






www.portable-ultrasound-scanner

Quality, Professional, Innovational


Specification

Master chip	dedicated master DSP
Magnification	5X ~ 200X
Take pictures / videos	built-in
Auxiliary light source	8 white LED lights
Static resolution	640x480
Imaging distance	manual adjustment 0 ~ 40mm
Fixed base	special base
Software	Professional intelligent skin analysis and testing software
Support system	Windows 7 64-bit, Windows 8 64-bit, Windows 10 64-bit
Connection	USB interface
Power supply	USB (5V DC)
Dynamic Frames	30f / s Under 600 LUX Brightness
Illumination range	0 ~ 30000LUX wire control adjustable
Supported languages	Simplified Chinese, Traditional Chinese, English
Product color	pearl white, other colors can be customized
Main body size	144mm (L) * 37mm (D)
Standard parts list	<ol style="list-style-type: none"> 1. USB high-definition skin analysis detector 2. Dedicated base 3. U key 4. Software 5. Paper manual 6. Black outer packing

Software interface



Images



Wuxi Biomedical Technology Co., Ltd. www.portable-ultrasound-scanner.com

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cnbiomedical.com

Mobile phone: 86-13382885910

BIE is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.